PROCEDIMENTO

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BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFÍA/BIBLIOGRAFIA

1. Sacks, D.B., *Carbohydrates*. *Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A, Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus*. *Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycemia*, (2006).
5. Trinder, P., *Ann. Clin. Biochem.*, (1969), 6, 24.
6. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/ LAB/99.1 Rev.2
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCC Press, (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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English - EN

INTENDED USE

GLUCOSE PAP SL is an *in vitro* diagnostic reagent intended for the quantitative determination of glucose in human serum and plasma samples on analyzers or semi-automatic analyzers. The standard is intended for the calibration of the reagent.
These *in vitro* diagnostic devices are for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Glucose is the main source of energy for the human body. Glucose is converted either into glycogen or into triglycerides to be stored. Glucose blood level is mainly regulated by two antagonist hormones: insulin and glucagon.
Glycemia disorders appear mostly in type I or type II diabetes as well as in gestational diabetes.

They can be also associated to various endocrinal, pancreatic or hepatic disorders, or linked to drugs. In normal health condition, glucose is filtered then reabsorbed by kidneys and is therefore not present in urine. Elevated concentrations in urine are observed when the blood concentration is high or in case of impaired tubular reabsorption. Glucose measurement in blood is indicated for diabetes, in screening, diagnosis or follow-up of patients. It is also indicated to monitor patients with symptoms of hyperglycemia or hypoglycemia.


LIMITATION OF USE

For diabetes assessment, the collection conditions and interpretation of serum glucose concentrations should follow local recommendations such as those published by WHO.⁽⁴⁾

The quantitative assay of glucose alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁵⁾

Enzymatic / PAP - End Point.

COMPOSITION

Reagent: R

Phosphate buffer, pH 7.4

Phenol	10	mmol/L
4-Aminoantipyrine	0.3	mmol/L
Glucose Oxidase	≥ 10 000	U/L
Peroxidase	≥ 700	U/L
Sodium azide	< 0.1	% (w/w)
Standard: Std (ref : GPSL-0497/0507/0707)		
D-Glucose	100	mg/dL
	5.55	mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- Reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

STABILITY

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels. The standard should be immediately and tightly capped to prevent contamination and evaporation. **On board stability :** The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES

Specimen ⁽¹⁾

- Serum.
- Plasma (lithium heparin).
- Plasma (Sodium fluoride / potassium oxalate (glycolysis inhibitors)).
- ☛ - CSF (Cerebrospinal Fluid) samples are compatible with this device based on data obtained by VitalScientific. It is recommended the use of CSF be validated by the laboratory.
- Using any other specimen type must be validated by the laboratory prior to use.
- Warnings and precautions**
- In samples collected without glycolysis inhibitors, blood cells must be removed rapidly to prevent the glucose loss (5% to 7% per hour in whole blood at room temperature).⁽¹⁾

- PAP method is not suitable for glucose measurement in urine due to the large amount of interfering endogenous substances present in this matrix. ^(1,2)
- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.
- Storage and stability ^(1,6)**
- Serum / Plasma (lithium heparin)
- 8 hours at room temperature
- 3 days at 2-8°C
- Plasma (sodium fluoride / potassium oxalate)
- 2 days at room temperature
- 7 days at 2-8°C

REFERENCE VALUES ⁽⁶⁾

Serum/plasma	mg/dL	mmol/L
Neonates	30 – 60	1.7 – 3.3
Children	60 – 100	3.3 – 5.6
Adults 18-60 y/o	74 – 106	4.1 – 5.9
Adults 60-90 y/o	82 – 115	4.6 – 6.4

***Note :** The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

☛INSTALLATION AND USE

For use on Selectra Pro analyzers :

- Consult operator manual.

- Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer. Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

PROCEDURE

Manual Procedure

<i>Wavelength :</i>	505 nm
<i>Optical path :</i>	1 cm
<i>Sample/ Reagent ratio :</i>	1:100
<i>Temperature :</i>	37 °C
<i>Read against reagent blank.</i>	

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 10 minutes.

Automatic Procedure

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

A Sample	x n	n = Calibrator/ standard
A Calibrator/ Standard		concentration

Conversion factor : mg/dL x 0.0555 = mmol/L

☛CALIBRATION

ELICAL 2 and Glucose Standard 100 mg/dL are traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay. Controls have to be performed :

- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

20.0 - 400.0 mg/dL (1.11 - 22.20 mmol/L) Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed.

This procedure extends the measuring range up to 2 000.0 mg/dL (111.01 mmol/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

LoD = 0.2 mg/dL (0.01 mmol/L)
LoQ = 10.0 mg/dL (0.56 mmol/L)

- Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean	Within-run	Total
	n	mg/dL	mmol/L	CV (%)
Level 1	80	37.4	2.08	0.7
Level 2	80	113.1	6.28	0.5
Level 3	80	284.0	15.76	0.7

- Correlation

A comparative study has been performed between GLUCOSE PAP SL reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples. The sample concentrations ranged from 22.2 to 384.9 mg/dL (1.23 - 21.36 mmol/L). The results are as follows :
Correlation coefficient : (r) = 1,000
Linear regression: y = 0.989x + 1.1 mg/dL (0.06 mmol/L)

- Limitations/Interferences

- Studies have been performed to determine the level of interference from different compounds. The following glucose levels were tested: 36.0 mg/dL, 108.1 mg/dL and 400.0 mg/dL. No significant interference is defined by a recovery ±10% of the initial value.
Unconjugated bilirubin: No significant interference up to 6.0 mg/dL (103 µmol/L).
Conjugated bilirubin: No significant interference up to 5.9 mg/dL (101 µmol/L).
Hemoglobin: No significant interference up to 300 mg/dL.
Triglycerides: No significant interference up to 920 mg/dL (10.40 mmol/L).
Ascorbic acid: No significant interference up to 2.0 mg/dL. Concentrations above the therapeutic levels may interfere and cause erroneous results.
Uric acid: No significant interference up to 23.0 mg/dL (1 368 µmol/L).
Methyl_dopa: No significant interference up to 0.8 mg/dL.
L-Dopa: Induces falsely low results at therapeutic concentrations.
Tolazamide: No significant interference up to 40.0 mg/dL.
Acetaminophen: No significant interference up to 30.0 mg/dL.
Do not use icteric or hemolyzed samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁷⁾

- Results can be falsely lowered by significant levels in the sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metamlzole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽⁸⁻⁹⁾

- On board stability/Calibration frequency

On Board Stability: 28 days

Calibration frequency: 28 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.
☛ *The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.*

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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