

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium hydroxide - reagent grade, ≥98%, pellets (anhydrous)

 Product Number:
 \$5881

 CAS Number:
 1310-73-2

 MDL:
 MFCD00003548

Formula Weight: 40.00

Formula Weight: 40.00 g/mol

NaOH

TEST Specification

Appearance (Color) White

Appearance (Form) Pellets
Solubility (Color) Colorless
Solubility (Turbidity) Clear

100 mg/ml, H2O

Titration with HCl \geq 98 % Impurity < 1.0 %

As Na2CO3

Recommended Retest Period -------

2 years

Specification: PRD.1.ZQ5.10000007758

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.