

# Flexibility for Your Vessels – Eppendorf SmartBlock™



SmartBlock removal from Eppendorf ThermoMixer® C



**condens.protect®**

Compatible Eppendorf SmartBlocks feature the condens.protect® symbol.

## Using different vessel sizes within your lab workflow?

Flexibility has never been this easy. Eppendorf offers a broad range of SmartBlocks for tubes from 0.2 mL to 50 mL as well as plates (MTP, DWP, and PCR plates 96/384). All SmartBlocks are equipped with the outstanding Eppendorf QuickRelease system that makes the block exchange super fast and easy. Just press the lever on the front of the block and the block can be removed – no tools are needed, done in seconds.

## Benefits

- > Optimized block design supports optimal temperature transfer within your sample to receive excellent results
- > Eppendorf SmartBlock-specific calibration by integrated sensor technology maximizes temperature accuracy and homogeneity
- > Every SmartBlock is individually inspected in production for temperature accuracy and is delivered incl. a serial-number specific certificate for documentation purpose
- > Insulated SmartBlocks for ergonomic operation – you won't burn your fingers
- > Automatic block recognition for safe handling

## Connect your ThermoMixer to VisioNize® Lab Suite



### Real-Time Monitoring

Monitor all device-related parameters like temperature and agitation, device interactions like set-point changes, and check the current status of your experiment or if a device is available



### Instant Notifications

Receive alarm and event notifications via email and SMS and always stay informed



### Easy Data Access

Export parameter, event data, and any user interactions

"Visionize Lab Suite helps us to automatically log device data and use it for further purposes. This makes our evaluations extremely robust and our prototyping complete."

Process Development Drug Substance,  
CureVac AG, Tübingen, Germany



**VisioNize** Lab Suite

**Sign up for free!**

Secure your samples with VisioNize Lab Suite for free!  
Contact us for more information: [www.eppendorf.com/visionize](http://www.eppendorf.com/visionize)

