

Specification Sheet

Product Name
Lactophenol blue solution for microscopy, for staining molds

Product Number

ZA.002.885

Product Brand

SIAL

CAS Number

Molecular Weight

TEST	SPECIFICATION
------	---------------

APPEARANCE (COLOR)	Blue to Very Dark Blue
--------------------	------------------------

APPEARANCE (FORM)	Liquid
-------------------	--------

DENSITY D20/4	1.145 - 1.174
---------------	---------------

WAVELENGTH (1) (UV)	601 - 611 nm
---------------------	--------------

EXTINCTION (1%/1CM) (1) (UV)	--
------------------------------	----

TEST	SPECIFICATION
------	---------------

WAVELENGTH (2) (UV)	313 - 323 nm
---------------------	--------------

EXTINCTION (1%/1CM) (2) (UV)	--
------------------------------	----

SOLVENT (UV)	ETOH
--------------	------

RECOMMENDED RETEST PERIOD	30 MONTHS
---------------------------	-----------